

TECHNICAL DATA SHEET

Recombinant Mouse CCL4 (MIP-1 β) (Carrier-free)

Catalog Number: 21-7098

RPx-Pro™ Recombinant Protein

PRODUCT INFORMATION

CONTENTS

Recombinant Mouse CCL4 (MIP-1 β) (Carrier-free)

DESCRIPTION

CCL4 (MIP-1 β) is a mitogen-inducible monokine and member of the CC chemokine family. It is a secreted protein that plays a role in the acute inflammatory response and is a chemoattractant for monocytes, neutrophils, eosinophils, NK cells, and CD4+ T lymphocytes. MIP-1 beta can signal through CCR5 and, by inducing downregulation of this receptor, can block HIV-1 infection of T cells.

MOLECULAR MASS

Recombinant Mouse MIP-1 beta consists of 69 amino acid residues, including the four highly conserved cysteine residues present in CC chemokines, and has a molecular weight of 7.8 kDa.

AMINO ACID SEQUENCE

APMGSDPPTS CCFSYTSRQL HRSFVMDYYE TSSLCSPAV VFLTKRGRQI CANPSEPWVT EYMSDLELN

SOURCE

E. coli

APPLICATIONS

Bioassay

PURITY

95 %

STORAGE

-20°C

PROTEIN CONTENT

Verified by UV Spectroscopy and/or SDS-PAGE gel.

ENDOTOXIN LEVEL

Endotoxin level is <0.1 ng/ μ g of protein (<1 EU/ μ g).

AUTHENTICITY

Verified by N-terminal and Mass Spectrometry analyses (when applicable).

CROSS REACTIVITY

Human

BIOACTIVITY

Using a concentration range of 5.0-20.0 ng/ml, activity is determined by its ability to chemoattract human blood monocytes.

RESEARCH AREAS

AIDS/HIV, Chemotaxis, Immune System, Inflammation, Neurobiology, Wound Healing

RECONSTITUTION

See Certificate of Analysis (COA) for lot specific reconstitution information.

REFERENCES

Schall TJ, Bacon K, Camp RD, Kaspari JW and Goeddel DV. 1993. J Exp Med. 177(6): 1821-1826. Menten P1, Wuyts A and Van Damme J. 2002. Cytokine Growth Factor Rev. 13(6): 455-481. Mueller A and Strange PG. 2004. Int J Biochem Cell Biol. 36(1): 35-38.

Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.