

TECHNICAL DATA SHEET

Recombinant Mouse PDGF-BB (Carrier-free)

Catalog Number: 21-7120

RPx-Pro™ Recombinant Protein

PRODUCT INFORMATION

CONTENTS

Recombinant Mouse PDGF-BB (Carrier-free)

DESCRIPTION

Platelet-derived growth factor-BB (PDGF-BB) is one of the five dimeric proteins belonging to the PDGF family. Other members are homodimers PDGF-AA, PDGF-CC, PDGF-DD, and the heterodimer PDGF-AB. The PDGFs are produced by platelets, stored in platelet alpha-granules and are released upon platelet activation. PDGF proteins are generally considered as potent mitogens for connective tissue cells, and can also be chemotactic for fibroblasts, smooth muscle cells, neutrophils and mononuclear cells. PDGF-BB plays a role in angiogenesis and smooth muscle formation during embryonic development, as well as wound healing and regeneration in adults. Two receptors have been identified for the PDGF proteins - PDGFR alpha and PDGFR beta. PDGF-BB binds to PDGFR alpha, as well as PDGFR beta which appears to be specific for PDGF-BB and PDGF-AB.

MOLECULAR MASS

Recombinant Murine PDGF-BB is a 24.4 kDa disulfide-linked homodimer made up of two B chains, each consisting of 109 amino acids.

AMINO ACID SEQUENCE

SLGSLAAAEPAVIAECKTRTEVFQISRNLIDRTNANFLVWPPCVEVQRCSGCCNNRNVCQRASQVQMRPVQVRKIEIVRKKPIFKKATVTL
LEDHLACKCE TIVTPRPVT

SOURCE

E. coli

APPLICATIONS

Bioassay

PURITY

98 %

STORAGE

-20°C

PROTEIN CONTENT

Verified by UV Spectroscopy and/or SDS-PAGE gel.

ENDOTOXIN LEVEL

Endotoxin level is <0.1 ng/μg of protein (<1 EU/μg).

AUTHENTICITY

Verified by N-terminal and Mass Spectrometry analyses (when applicable).

CROSS REACTIVITY

Monkey, Human

BIOACTIVITY

The expected ED₅₀ is 1.0-2.0 ng/ml, as determined by the dose-dependent stimulation of BALB/c 3T3 cell proliferation.

RESEARCH AREAS

Angiogenesis & Cardiovascular, Bone and Cartilage, Cancer, Differentiation, Immune System, Neurobiology, Stem Cells, Wound Healing

RECONSTITUTION

See Certificate of Analysis (COA) for lot specific reconstitution information.

REFERENCES

Fredriksson L, Li H and Eriksson U. 2004. Cytokine Growth Factor Rev. 15(4): 197-204. Heldin CH and Westermark B. 1999. Physiol Rev. 79(4): 1283-1316. Kamimura M, Bea F, Akizawa T, Katus HA, Kreuzer J and Viedt C. 2004. Hypertension. 44(6): 944-951. Anand-Apte B and Zetter B. 1997. Stem Cells. 15(4): 259-267. Siegbahn A, Hammacher A, Westermark B and Heldin CH. 1990. J Clin Invest. 85(3): 916-920. Andrae J, Gallini R and Betsholtz C. 2008. Genes Dev. 22(10): 1276-1312.

Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.