

TECHNICAL DATA SHEET

Recombinant Human Adiponectin (Acrp30) (Carrier-free)

Catalog Number: 21-7165

RPx-Pro™ Recombinant Protein

PRODUCT INFORMATION

CONTENTS

Recombinant Human Adiponectin (Acrp30) (Carrier-free)

DESCRIPTION

Adiponectin, or Acrp30, is a secreted hormone protein (adipokine) that is expressed exclusively by adipocytes and is present in circulating plasma. Low plasma levels of this protein are associated with insulin resistance and obesity. Adiponectin signals through the AdipoR1 and AdipoR2 receptors, and its interaction with various growth factors may, in part, account for a role in the regulation of cell growth and tissue remodeling. This protein had been shown to exhibit both pro- and anti-inflammatory effects.

MOLECULAR MASS

Recombinant Human Adiponectin, a glycosylated multimer, has a calculated molecular weight of 25.9 kDa. It contains amino acids Glu-19 to Asn-244 of the precursor protein fused to an N-terminal histidine tag. Under reducing conditions, the glycosylated monomer migrates at about 35 kDa by SDS PAGE.

AMINO ACID SEQUENCE

RGHHHHHHHH ETTTQGPGLV LPLPKGACTG WMAGIPGHPG HNGAPGRDGR DGTGPEKGEK GDPGLIGPKG DIGETGVPGA EGPRGFPGIQ
GRKGEPGEGA YVYRSAFSVG LETYVTIPNM PIRFTKIFYN QQNHYDGSTG KFHCNIPGLY YFAYHITVYM KDVK/SLFKK DKAMLFTYDQ YQENNVDAQS
GSVLLHLEVG DQVWLQVYGE GERNGLYADN DNDSTFTGFL LYHDTN

SOURCE

Hi-5 Insect cells

APPLICATIONS

Bioassay

PURITY

97 %

STORAGE

-20°C

PROTEIN CONTENT

Verified by UV Spectroscopy and/or SDS-PAGE gel.

ENDOTOXIN LEVEL

Endotoxin level is <0.1 ng/µg of protein (<1 EU/µg).

AUTHENTICITY

Verified by N-terminal and Mass Spectrometry analyses (when applicable).

CROSS REACTIVITY

Mouse, Rat

BIOACTIVITY

The expected ED₅₀ is 3.0-6.0 µg/ml, as determined by a cytotoxicity assay using murine M1 cells.

RESEARCH AREAS

Diabetes, Inflammation

RECONSTITUTION

See Certificate of Analysis (COA) for lot specific reconstitution information.

REFERENCES

Nedvídková J, Smitka K, Kopský V and Hainer V. 2005. *Physiol Res.* 54(2): 133-140. Díez JJ and Iglesias P. 2003. *Eur J Endocrinol.* 148(3): 293-300. Wang Y, Lam KS, Xu JY, Lu G, Xu LY, Cooper GJ and Xu A. 2005. *J Biol Chem.* 280(18): 18341-18347. Yamauchi T, Nio Y, Maki T, Kobayashi M, Takazawa T, Iwabu M, Okada-Iwabu M, Kawamoto S, Kubota N, Kubota T, Ito Y, Kamon J et al. 2007. *Nat Med.* 13(3): 332-339.

Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.