

## TECHNICAL DATA SHEET

# Recombinant Human EG-VEGF (Carrier-free)

Catalog Number: 21-7179

## RPx-Pro™ Recombinant Protein

### PRODUCT INFORMATION

#### CONTENTS

Recombinant Human EG-VEGF (Carrier-free)

#### DESCRIPTION

EG-VEGF, also known as prokineticin-1, is a member of the AVIT (prokineticin) family known to strongly contract gastrointestinal (GI) smooth muscle. EG-VEGF is a secreted angiogenetic mitogen growth factor that induces proliferation, migration, and fenestration in capillary endothelial cells derived from endocrine glands. The steroidogenic glands (ovary, testis, adrenal gland, and placenta) express EG-VEGF and expression is known to be induced by hypoxia. EG-VEGF and VEGF may function in a coordinated manner as their expression is often complementary. EG-VEGF is associated with certain cancers and is known to promote neuroblastoma progression. The human EG-VEGF gene codes for a 105 amino acid polypeptide containing an N-terminal signal sequence of 19 amino acids.

#### MOLECULAR MASS

Recombinant Human EG-VEGF is a 9.6 kDa protein consisting of 86 amino acid residues, including ten cysteine residues that potentially form five pairs of intra-molecular disulfide bonds.

#### AMINO ACID SEQUENCE

AVITGACERD VQCGAGTCCA ISLWLRGLRM CTPLGREGEE CHPGSHKVPF FRKRKHTTCP CLPNLLCSRF PDGRYRCSMD LKNINF

#### SOURCE

E.coli

#### APPLICATIONS

Bioassay

#### PURITY

98 %

#### STORAGE

-20°C

#### PROTEIN CONTENT

Content Verified by UV Spectroscopy and/or SDS-PAGE gel.

#### ENDOTOXIN LEVEL

Endotoxin level is <0.1 ng/μg of protein (<1EU/μg).

#### AUTHENTICITY

Verified by N-terminal and Mass Spectrometry analyses (when applicable).

#### CROSS REACTIVITY

Human

#### BIOACTIVITY

Data not available at this time.

#### RESEARCH AREAS

Angiogenesis & Cardiovascular, Proliferation

#### RECONSTITUTION

See Certificate of Analysis (COA) for lot specific reconstitution information.

#### REFERENCES

Cook IH, Evans J, Maldonado-Pérez D, Critchley HO, Sales KJ and Jabbour HN. 2010. Mol Hum Reprod. 16(3):158-69. Ngan ES, Sit FY, Lee K, Miao X, Yuan Z, Wang W, Nicholls JM, Wong KK, Garcia-Barcelo M, Lui VC and Tam PK. 2007. Clin Cancer. Res. 13(3):868-75. Denison FC, Battersby S, King AE, Szuber M and Jabbour HN. 2008. Endocrinology. 149(7):3470-7. Waddell JM, Evans J, Jabbour HN and Denison FC. 2011. Hum Reprod. 26(1):67-75. Gorowiec MR, Catalano RD, Norman JE, Denison FC and Jabbour HN. 2011. Am J Pathol. 179(6):2709-19.

Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.