

TECHNICAL DATA SHEET

# Recombinant Mouse CCL2 (MCP-1)

Catalog Number: 21-8152

## RPx-Pro™ Recombinant Protein PRODUCT INFORMATION

### CONTENTS

Recombinant Mouse CCL2 (MCP-1)

### DESCRIPTION

Monocyte chemoattractant protein-1, MCP-1, also known as Chemokine (C-C motif) ligand 2, CCL2, is a member of the CC chemokine family and is expressed in endothelial cells, smooth muscle cells, and monocytes. At sites of inflammation MCP-1 is known to recruit monocytes, memory T cells and dendritic cells. MCP-1 is implicated in the pathogenesis of psoriasis, rheumatoid arthritis and atherosclerosis.

### MOLECULAR MASS

Recombinant Mouse MCP-1 (CCL2) is a 13.8 kDa protein containing 125 amino acid residues

### AMINO ACID SEQUENCE

QPDAVNAPLT CCYSFTSKMI PMSRLESYKR ITSSRCPKEA VVFTVTKLKRE VCADPKKEWV QTYIKNLDRN QMRSEPTTLF KTAS-  
ALRSSA PLNVKLTRKS EANASTTFST TTSSTSVGVT SVTVN

### SOURCE

E. Coli

### PROTEIN CONTENT

Content Verified by UV Spectroscopy and/or SDS-PAGE gel.

### PURITY

98 %

### ENDOTOXIN LEVEL

Endotoxin level is <0.1 ng/μg of protein (<1EU/μg).

### AUTHENTICITY

Verified by N-terminal and Mass Spectrometry analyses (when applicable).

### CROSS REACTIVITY

Bacteria, Human, Mouse, Rabbit, Monkey

### BIOACTIVITY

Determined by its ability to chemoattract Balb/C mouse spleen MNCs using a concentration range of 1.0-20.0 ng/ml.

### APPLICATIONS

Bioassay

### RESEARCH AREAS

Angiogenesis/Cardiovascular, Cancer, Chemotaxis, Immune System, Inflammation, Neurobiology, Wound Healing, Allergy.

### STORAGE

-20°C

### RECONSTITUTION

See Certificate of Analysis (COA) for lot specific reconstitution information

### REFERENCES

Carr, M. W.; Roth, S. J.; Luther, E.; Rose, S. S.; Springer, T. A. (1994). "Monocyte chemoattractant protein 1 acts as a T-lymphocyte chemoattractant". *Proceedings of the National Academy of Sciences of the United States of America* 91 (9): 3652–3656. Xu, L. L.; Warren, M. K.; Rose, W. L.; Gong, W.; Wang, J. M. (1996). "Human recombinant monocyte chemoattractant protein and other C-C chemokines bind and induce directional migration of dendritic cells in vitro". *Journal of leukocyte biology* 60 (3): 365–371. Xia, M.; Sui, Z. (2009). "Recent developments in CCR2 antagonists". *Expert Opinion on Therapeutic Patents* 19 (3): 295–303.

NOTE: Please choose the appropriate format for each application. Citations are provided as a convenience to you; please consult Materials and Methods sections for additional details about the use of any product in these publications.

### For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.