

TECHNICAL DATA SHEET

Recombinant Human CXCL9 (MIG) (Carrier-free)

Catalog Number: 21-8299

RPx-Pro™ Recombinant Protein

PRODUCT INFORMATION

CONTENTS

Recombinant Human CXCL9 (MIG) (Carrier-free)

DESCRIPTION

Monokine induced by gamma interferon, MIG, also known as chemokine (C-X-C motif) ligand 9, CXCL9. MIG is a T-cell chemoattractant belonging to the CXC chemokine family inducing effects through interaction with chemokine receptor CXCR3. Through interaction with activated CD4+ Th1 cells, CD8+ T cells, IL-2 activated T lymphocytes, and NK cells, MIG plays an important role in leukocyte trafficking.

MOLECULAR MASS

Recombinant Human MIG (CXCL9) is a 11.7 kDa protein composed of 103 amino acid residues.

AMINO ACID SEQUENCE

TPVVRKGRCS CISTNQGTIH LQSLKDLKQF APSPSCEKIE IIATLKNQVQ TCLNPDSADV KELIKKWEKQ VSQKKKQKNG KKHQKKKVLK VRKSQRSRQK KTT

SOURCE

E. Coli

APPLICATIONS

Bioassay

PURITY

98 %

STORAGE

-20°C

PROTEIN CONTENT

Content Verified by UV Spectroscopy and/or SDS-PAGE

ENDOTOXIN LEVEL

Endotoxin level is <0.1 ng/μg of protein (<1 EU/μg).

AUTHENTICITY

Verified by N-terminal and Mass Spectrometry analyses (when applicable).

CROSS REACTIVITY

N/A

BIOACTIVITY

Determined by its ability to chemoattract human peripheral blood T lymphocytes using a concentration range of 10.0-100.0 ng/ml.

RESEARCH AREAS

Angiogenesis/Cardiovascular; Chemotaxis; Immune System; Inflammation ; Wound Healing; Transplantation

RECONSTITUTION

See Certificate of Analysis (COA) for lot specific reconstitution information.

REFERENCES

Tensen CP, Flier J, Van Der Raaij-Helmer EM, Sampat-Sardjoepersad S, Van Der Schors RC, Leurs R, Scheper RJ, Boorsma DM and Willemze R 1999 J Invest Dermatol 112(5): 716–22. Weng Y, SJ Siciliano, KE Waldburger, A Sirotna-Meisher, MJ Staruch, BL Daugherty, SL Gould, MS Springer and JA DeMartino 1998 J Biol Chem 273: 18288–18291.

Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.