

TECHNICAL DATA SHEET

Recombinant Mouse TNF- α (Carrier-free)

Catalog Number: 21-8321

RPx-Pro™ Recombinant Protein
PRODUCT INFORMATION

CONTENTS

Recombinant Mouse TNF- α (Carrier-free)

DESCRIPTION

TNF-alpha is a clinically important cytokine with strong pro-inflammatory activities. TNF-alpha is secreted by a variety of immune and non-immune cells including adipocytes, activated monocytes, macrophages, B cells, T cells and fibroblasts. TNF-alpha signals through two receptors, TNFR1 and TNFR2. TNF-alpha has strong cytolytic activities against a wide variety of tumor cells and is critically important in combat against bacterial infections. The inflammatory effects of TNF-alpha are also deleterious in causation of septic shock and autoimmune diseases including rheumatoid arthritis, inflammation, and diabetes.

MOLECULAR MASS

Recombinant Mouse TNF- α is a 17.3 kDa 156 amino acid protein representing the full length transmembrane protein.

AMINO ACID SEQUENCE

MLRSSSQNSS DKPVAHVVAN HQVEEQLEWL SQRANALLAN GMDLKDNLV VPADGLYLVY SQVLFGQGC PDYVLLTHTV
SRFAISYQEK VNLLSAVKSP CPKDTPEGAE LKPWYEPIYL GGVFQLEKGD QLSAEVNLPK YLDFAESGQV YFGVIAL

SOURCE

E. Coli

APPLICATIONS

Bioassay

PURITY

98 %

STORAGE

-20°C

PROTEIN CONTENT

Content Verified by UV Spectroscopy and/or SDS-PAGE

ENDOTOXIN LEVEL

Endotoxin level is <0.1 ng/ μ g of protein (<1 EU/ μ g).

AUTHENTICITY

Verified by N-terminal and Mass Spectrometry analyses (when applicable).

CROSS REACTIVITY

Human, Monkey, Rat, Cow

BIOACTIVITY

ED₅₀ is \leq 0.1 ng/ml, corresponding to a specific activity of $\geq 71 \times 10^7$ units/mg calculated by the cytolysis of mouse L929 cells in the presence of actinomycin D.

RESEARCH AREAS

AIDS/HIV; Angiogenesis/Cardiovascular; Apoptosis; Cancer; Diabetes / Weight Regulation; Immune System; Inflammation; Neurobiology; Stem Cells & Differentiation; TNF Superfamily; Allergy; Transplantation

RECONSTITUTION

See Certificate of Analysis (COA) for lot specific reconstitution information.

REFERENCES

Feldmann M and RN Maini 2010 J Immunol 185: 791-794. Palladino MA, FR Bahjat, EA Theodorakis and LL Moldawer 2003 Nat Rev Drug Discov 2: 736-46. Swardfager W, K Lancôt, L Rothenburg, A Wong, J Cappel and N Herrmann 2010 Biol Psychiatry 68: 930-941. Locksley RM, N Killeen and MJ Lenardo 2001 Cell 104: 487-501. Kolb WP and GA Granger Proc Natl Acad Sci USA 61(4): 1250-1255. Ruddle NH and BH Waksman 1968 J Exp Med 128(6): 1267-1279. Carswell EA, LJ Old, RL Kassel, S Green, N Fiore and B Williamson 1975 Proc Natl Acad Sci USA 72: 3666-3670.

Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.