

**TECHNICAL DATA SHEET**

**Recombinant Mouse LIX (70 a.a.) (CXCL5) (Carrier-Free)**

Catalog Number: 21-9141

**RPx-Pro™ Recombinant Protein**  
**PRODUCT INFORMATION**

**CONTENTS**

Recombinant Mouse LIX (70 a.a.) (CXCL5) (Carrier-Free)

**DESCRIPTION**

CXCL6, also known as LIX in mice, is a CXC chemokine that signals through the CXCR2 receptor. It is expressed in monocytes, platelets, endothelial cells, and mast cells. LIX is a chemoattractant for neutrophils. The two naturally occurring variants of LIX, LIX 1-78 (GCP-2) and LIX 9-78 (GCP-2), contain 78 and 70 amino acid residues, respectively. LIX contains the four conserved cysteine residues present in CXC chemokines, and also contains the 'ELR' motif common to CXC chemokines that binds to the CXCR1 and CXCR2 receptors.

**MOLECULAR MASS**

Recombinant Mouse LIX is a 7.6 kDa protein consisting of amino acid residues 9-78 of the full length LIX protein (70 amino acids).

**AMINO ACID SEQUENCE**

TELRVCVCLTV TPKINPKLIA NLEVIPAGPQ CPTVEVIAKL KNQKEVCLDP EAPVIKKIIQ KILGSDKKKA

**SOURCE**

E. coli

**APPLICATIONS**

Bioassay

**PURITY**

98 %

**STORAGE**

-20°C

**PROTEIN CONTENT**

Content Verified by UV Spectroscopy and/or SDS-PAGE gel.

**ENDOTOXIN LEVEL**

Endotoxin level is <0.1 ng/µg of protein (<1EU/µg).

**AUTHENTICITY**

Verified by N-terminal and Mass Spectrometry analyses (when applicable).

**CROSS REACTIVITY**

Human, Mouse

**BIOACTIVITY**

Determined by its ability to chemoattract human neutrophils using a concentration range of 10.0-100.0 ng/ml.

**RESEARCH AREAS**

Angiogenesis/Cardiovascular, Wound Healing, Chemotaxis, Inflammation

**RECONSTITUTION**

See Certificate of Analysis (COA) for lot specific reconstitution information.

**REFERENCES**

Mei, J. CXCL5 regulates chemokine scavenging and pulmonary host defense to bacterial infection. 2010. *Immunity*; 33(1):106-17. Hosking, M.P. A protective role for ELR+ chemokines during acute viral encephalomyelitis. 2009. *PLoS Pathogens*; 5(11):e1000648.

Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.