

TECHNICAL DATA SHEET

Recombinant Human BD-4 (Carrier-Free)

Catalog Number: 21-9174

RPx-Pro™ Recombinant Protein
PRODUCT INFORMATION

CONTENTS

Recombinant Human BD-4 (Carrier-Free)

DESCRIPTION

Defensins (alpha and beta) are cationic peptides with a broad spectrum of antimicrobial activity that comprise an important arm of the innate immune system. The α -defensins are distinguished from the β -defensins by the pairing of their three disulfide bonds. To date, six human β -defensins have been identified; BD-1, BD-2, BD-3, BD-4, BD-5 and BD-6. β -defensins are expressed on some leukocytes and at epithelial surfaces. In addition to their direct antimicrobial activities, they can act as chemoattractants towards immature dendritic cells and memory T cells. The β -defensin proteins are expressed as the C-terminal portion of precursors, and are released by proteolytic cleavage of a signal sequence and, in some cases, a propeptide sequence. β -defensins contain a six-cysteine motif that forms three intra-molecular disulfide bonds. BD-4 is expressed in the testes, stomach, uterus, neutrophils, thyroid, lungs and kidneys. In addition to its direct antimicrobial activities, BD-4 is chemoattractant towards human blood monocytes.

MOLECULAR MASS

Recombinant Human BD-4 is a 6.0 kDa protein containing 50 amino acid residues.

AMINO ACID SEQUENCE

EFELDRICGY GTARCRKKCR SQEYRIGRCP NTYACCLRKW DESLLNRTKP

SOURCE

E. coli

APPLICATIONS

Bioassay

PURITY

98 %

STORAGE

-20°C

PROTEIN CONTENT

Content Verified by UV Spectroscopy and/or SDS-PAGE gel.

ENDOTOXIN LEVEL

Endotoxin level is <0.1 ng/μg of protein (<1EU/μg).

AUTHENTICITY

Verified by N-terminal and Mass Spectrometry analyses (when applicable).

CROSS REACTIVITY

N/A

BIOACTIVITY

Determined by its ability to chemoattract human monocytes using a concentration range of 0.1-100.0 ng/ml.

RESEARCH AREAS

Immune System, Chemotaxis

RECONSTITUTION

See Certificate of Analysis (COA) for lot specific reconstitution information.

REFERENCES

Medina-Inojosa J, Somers V, Jenkins S, Zundel J, Johnson L, Grimes C, Lopez-Jimenez F. *Obes Open Access*. 2017;3(1). doi: 10.16966/2380-5528.127.
 Abdul-Malak O, Vodovotz Y, Zaaqoq A, Guardado J, Almahmoud K, Yin J, Zuckerbraun B, Peitzman AB, Sperry J, Billiar TR, Namas RA. *Mediators Inflamm*. 2016;2016:7950374.

Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.