

TECHNICAL DATA SHEET

# FITC Mouse IgG1 Isotype Control (MOPC-21)

Catalog Number: 35-4714

## PRODUCT INFORMATION

**Contents:** FITC Mouse IgG1 Isotype Control (MOPC-21)

**Isotype:** Mouse IgG1

**Concentration:** 0.5 mg/mL

**Clone:** MOPC-21

**Reactivity:** N/A

**Formulation:** 10 mM NaH<sub>2</sub>PO<sub>4</sub>, 150 mM NaCl, 0.09% NaN<sub>3</sub>,  
0.1% gelatin, pH7.2

## DESCRIPTION

The MOPC-21 immunoglobulin is useful as an isotype-matched control. The MOPC-21 immunoglobulin has an unknown binding specificity and is used as an isotype control for mouse IgG1 antibodies.

## PREPARATION & STORAGE

This monoclonal antibody preparation was purified from tissue culture supernatant via affinity chromatography. For In Vivo Ready™ (IVR) products, each preparation is also evaluated for endotoxin levels using the LAL assay. It is recommended to store the product undiluted at 4°C. Do not freeze.

## APPLICATION NOTES

This antibody preparation has been quality-tested for flow cytometry using an appropriate cell type. The isotype control should be used at the same concentration as the primary antibody for which it is a control.

NOTE: Please choose the appropriate format for each application. Citations are provided as a convenience to you; please consult Materials and Methods sections for additional details about the use of any product in these publications.

## For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biosciences Corporation. © 2013 Tonbo Biosciences.