

TECHNICAL DATA SHEET

Goat anti-Mouse IgG (H&L) - F(ab')₂ fragment, FITC Conjugate, Cross-Adsorbed

Catalog Number: 35-8057

PRODUCT INFORMATION

Contents: Goat anti-Mouse IgG (H&L) - F(ab')₂ fragment, FITC Conjugate, min x w/bovine, goat, human, rabbit, or rat IgG

Reactivity: Mouse

Specificity: Mouse IgG (H&L), F(ab')₂ Fragment

Host: Goat

Cross Adsorption: min x w/bovine, goat, human, rabbit, or rat IgG

Cross Reactivity: Based on IEP, no reactivity is observed to: · non-immunoglobulin mouse serum immunoglobulins · bovine, goat, human, rabbit or rat IgG

Concentration: 1.5 mg/mL

Formulation: 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free

Reconstitution:

Storage Conditions: Store at 2-8 °C

Use By: 1 year from date of receipt. Prepare working dilution prior to use and then discard.

Purity: > 90% based on SDS-PAGE small amounts of intact IgG may be present

Purification: Affinity purified using solid phase Mouse IgG

Country of Origin: Goat serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.

PREPARATION & STORAGE

This product was Affinity purified using solid phase Mouse IgG. Store at 2-8 °C

APPLICATION NOTES

This conjugate is suitable for flow cytometry, immunofluorescence, ELISA and immunomicroscopy. The optimal working dilution should be determined by the investigator. Suggested starting dilution(κ): · 1:20 - 1:2,000 for most applications Based on IEP, this antibody reacts with: · heavy (γ) chains on mouse IgG · light chains on all mouse immunoglobulins.

Tonbo Biosciences tests all antibodies by flow cytometry. Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.