

TECHNICAL DATA SHEET

In Vivo Ready™ Anti-Mouse CD3e (145-2C11)

Catalog Number: 40-0031

PRODUCT INFORMATION

Contents: In Vivo Ready™ Anti-Mouse CD3e (145-2C11)

Isotype: Armenian Hamster IgG

Concentration: 2 mg/mL

Clone: 145-2C11

Reactivity: Mouse

Formulation: 10 mM NaH₂PO₄, 150 mM NaCl, pH7.2

Endotoxin Level: Less than or equal to 0.01 EU/ug, as determined by the LaL assay

DESCRIPTION

The 145-2C11 antibody is specific for mouse CD3e, also known as CD3 epsilon, a 20 kDa subunit of the T cell receptor complex, along with CD3 gamma and CD3 delta. These integral membrane protein chains assemble with additional chains of the T cell receptor (TCR), as well as CD3 zeta chain, to form the T cell receptor – CD3 complex. Together with co-receptors CD4 or CD8, the complex serves to recognize antigens bound to MHC molecules on antigen-presenting cells. Such interactions promote T cell receptor signaling (T cell activation) and can result in a number of cellular responses including proliferation, differentiation, production of cytokines or activation-induced cell death. CD3 is differentially expressed during thymocyte-to-T cell development and on all mature T cells. The 145-2C11 antibody is a widely used phenotypic marker for mouse T cells. In addition, binding of 145-2C11 antibody to CD3e can induce cell activation. A recent publication of the crystal structure of a murine CD3e-mitogenic antibody complex provides further insight into the action of commonly used agonist antibodies (Fernandes, R.A. et al. 2012. J. Biol. Chem. 287: 13324-13335).

PREPARATION & STORAGE

This monoclonal antibody preparation was purified from tissue culture supernatant via affinity chromatography. For In Vivo Ready™ (IVR) products, each preparation is also evaluated for endotoxin levels using the LAL assay. It is recommended to store the product undiluted at 4°C. Do not freeze.

APPLICATION NOTES

This purified format is guaranteed to be >90% pure as determined by SDS-PAGE analysis. Citations are provided as a convenience to you - please consult Materials and Methods sections for additional details about the use of any product in these publications.

REFERENCES

Staepli R, Ludigs K, Heinz LX, Sequin-Estevez Q, Ferero I, Braun M, Schroder K, Rebsamen M, Tardivel A, Mattmann C, MacDonald HR, Romero P, Reith W, Guarda G, and Tschopp J. 2012. J. Immunol. 188: 3820-3828. (in vitro activation)Todo T, Wu G, Chai NN, He Y, Martins G, Gupta A, Fair J, Liu NY, Jordan S, and Klein A. 2012. Int. Immunol. 10:1093. (in vivo assay)Mira E, Leon B, Barber DF, Jimenez-Baranda S, Goya I, Almonacid L, Marquez G, Zaballo A, Martinez C, Stein JV, Ardavin C and Manes S. 2012. J. Immunol. (in vitro activation, immunohistochemistry – frozen tissue)Becker-Herman A, Meyer-Bahlburg A, Schwartz MA, Jackson SW, Hudkins KL, Liu C, Sather BD, Khim S, Liggitt D, Song W, Silverman GJ, Alpers CE and Rawlings DJ. 2011. J. Exp. Med. 208:2033-2042. (immunofluorescence microscopy – OCT embedded frozen tissue)Salmond RJ, Filby A, Pirinen N, Magee AI, and Zamoyska R. 2010. Blood. 117: 108-117. (immunoprecipitation)Tilley SL, Jaradat M, Stapleton C, Dixon D, Hua X, Erikson CJ, McCaskill JG, Chason KD, Liao G, Jania L, Koller BH, and Jetten AM. 2007. J. Immunol. 178: 3208-18. (immunohistochemistry – frozen tissue)Isakov N, Wange RL, Burgess WH, Watts JD, Aebersold R, and Samelson LE. 1995. J. Exp. Med. 181:375-380. (in vitro activation, immunoprecipitation)Salvadori S, Gansbacher B, Pizzimenti AM, and Zier KS. 1994. J. Immunol. 153: 5176 - 5182. (western blotting)Leo O, Foo M, Sachs D, Samuelson L, and Bluestone J. 1987. Proc. Natl. Acad. Sci. USA. 84: 1374 (Origination of clone 145-2C11, in vitro activation, in vitro blocking, immunoprecipitation)

NOTE: Please choose the appropriate format for each application. Citations are provided as a convenience to you; please consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.