

TECHNICAL DATA SHEET

In Vivo Ready™ Anti-Mouse CD40 (FGK45)

Catalog Number: 40-8050

PRODUCT INFORMATION

Contents: In Vivo Ready™ Anti-Mouse CD40 (FGK45)

Isotype: Rat IgG2a

Concentration: 2.0 mg/mL

Clone: FGK45

Reactivity: Mouse

Use By: 12 months from date of receipt

Storage Conditions: 2-8°C

Endotoxin Level: Less than or equal to 0.01 EU/ug, as determined by the LaL assay

Formulation: 10 mM NaH₂PO₄, 150 mM NaCl, pH7.2

DESCRIPTION

The FGK45 antibody reacts with mouse CD40, a 40-50 kDa glycoprotein member of the tumor necrosis factor receptor (TNFR) family. CD40 is constitutively expressed by antigen presenting cells, including dendritic cells, B cells, and macrophages. The interaction of CD40 and its ligand CD154, expressed on activated T cells, plays an essential role in mediating a broad range of immune and inflammatory responses including T cell dependent immunoglobulin class switching, germinal center formation, and the development of memory B cells.

The FGK45 clone is reported to be functionally useful for both in vitro and in vivo activation of CD40 expressing antigen presenting cells.

PREPARATION & STORAGE

This monoclonal antibody preparation was purified from tissue culture supernatant via affinity chromatography. For In Vivo Ready™ (IVR) products, each preparation is also evaluated for endotoxin levels using the LAL assay. It is recommended to store the product undiluted at 4°C. Do not freeze.

APPLICATION NOTES

Tonbo Biosciences tests all of our antibodies by flow cytometry. Citations may be provided as a resource for additional applications that have not been validated by Tonbo Biosciences - please consult Materials and Methods sections for additional details about the use of any product in these publications.

REFERENCES

Shepherd DM, Kerkvliet NI. 1999. J Immunol. Sep 1;163(5):2470-2477. Mauri C, Mars LT, Londei M. 2000. Nat Med. Jun;6(6):673-679. Hellings PW, Kasran A, Bullens D, Overbergh L, Mathieu C, Heremans H, Matthys P, Boon L, Jorissen M, Ceuppens JL. 2006. J Immunol. Oct 15;177(8):5138-5144. Bartholdy C, Kauffmann SO, Christensen JP, Thomsen AR. 2007. J Immunol. Feb 1;178(3):1662-1670.

Tonbo Biosciences tests all antibodies by flow cytometry. Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.