

TECHNICAL DATA SHEET

In Vivo Ready™ Anti-Mouse CD279 (PD-1) (J43.1)

Catalog Number: 40-9985

PRODUCT INFORMATION**Contents:** In Vivo Ready™ Anti-Mouse CD279**Isotype:** Armenian Hamster IgG**Concentration:** 2.0 mg/mL**Clone:** J43.1**Reactivity:** Mouse**Use By:** 12 months from date of receipt**Storage Conditions:** 2-8°C**Endotoxin Level:** Less than or equal to 0.01 EU/ug, as determined by the LaL assay**Formulation:** 10 mM NaH₂PO₄, 150 mM NaCl, pH7.2**DESCRIPTION**

The J43.1 antibody is specific for mouse CD279, also known as programmed death-1 (PD-1), a 55 kDa glycoprotein which can co-regulate T cell antigen receptor signaling and therefore modulate T cell activation. PD-1 exists in a monomeric form that is expressed by CD4- CD8- thymocytes, where it participates in the processes of clonal selection, elimination of autoreactive lymphocytes, and development of tolerance. PD-1 expression is also inducible upon activation of mature T cells, where it has been proposed to interact with the co-stimulatory receptor CD80 to limit T cell activation. Two ligands for PD-1, known as PD-L1 (B7-H1) and PD-L2 (B7-DC) are differentially expressed on T and B cells, monocytes, macrophages, NK cells or dendritic cells. PD-1 is a member of a family of receptors including CD28 and CTLA-4 (CD152), which interact with "B7" ligands to provide a balance of co-stimulatory /co-inhibitory signaling important in T cell activation, tolerance, and autoimmunity. The J43.1 antibody may be used as a marker for PD-1 expression, and is commonly used for analysis of receptor-ligand interaction and function(s) in vitro and in vivo.

PREPARATION & STORAGE

This monoclonal antibody preparation was purified from tissue culture supernatant via affinity chromatography. For In Vivo Ready™ (IVR) products, each preparation is also evaluated for endotoxin levels using the LAL assay. It is recommended to store the product undiluted at 4°C. Do not freeze.

APPLICATION NOTES

This purified format is guaranteed to be >90% pure as determined by SDS-PAGE analysis. Tonbo Biosciences tests all of our antibodies by flow cytometry. Citations may be provided as a resource for additional applications that have not been validated by Tonbo Biosciences - please consult Materials and Methods sections for additional details about the use of any product in these publications.

REFERENCES

Hams E, McCarron MJ, Amu S, Yagita H, Azuma M, Chen L, and Fallon PG. 2011. J. Immunol. 186:5648-5655. (in vivo blocking) Rivas MN, Weatherly K, Hazzan M, Vokaer B, Dremier S, Gaudray F, Goldman M, Salmon I, and Braun MY. 2009. 183:4284-4291. (in vitro blocking) Koehn BH, Ford ML, Ferrer IR, Borom K, Gangappa S, Kirk AD, and Larsen CP. 2008. J. Immunol. 181:5313-5322. (in vivo blocking) Brooks DG, Ha S-J, Elsaesser H, Sharpe AH, Freeman GJ, and Oldstone MBA. 2008. Proc. Natl. Acad. Sci. 105:20428-20433. (Flow cytometry) Ansari MJI, Salama AD, Chitnis T, Smith RN, Yagita H, Akiba H, Yamazaki T, Azuma M, Isai H, Khoury SJ, Auchincloss H, and Sayegh MH. 2003. J. Exp. Med. 198:63-71. (Immunohistochemistry – frozen tissue, in vivo blocking) Agata Y, Kawasaki A, Nishimura H, Ishida Y, Tsubat Y, Yagita H, and Honjo T. 1996. Int. Immunol. 8:765-772. (Immunoprecipitation)

Tonbo Biosciences tests all antibodies by flow cytometry. Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.