

TECHNICAL DATA SHEET

Purified Anti-Human CD40 (G28.5)

Catalog Number: 70-0410

PRODUCT INFORMATION

Contents: Purified Anti-Human CD40 (G28.5)

Isotype: Mouse IgG1, kappa

Concentration: 0.5 mg/mL

Clone: G28.5

Reactivity: Human

Use By: 12 months from date of receipt

Storage Conditions: 2-8°C

Formulation: 10 mM NaH₂PO₄, 150 mM NaCl, 0.09% NaN₃, pH7.2

DESCRIPTION

The G28.5 antibody reacts with human CD40, a 45-48 kDa glycoprotein that is a member of the TNF-receptor superfamily. CD40 is expressed on B cells, follicular dendritic cells, macrophages, endothelial cells, and fibroblasts. CD40 functions as a receptor for CD154 (CD40 ligand) and mediates a wide variety of immune and inflammatory responses including co-stimulation, B cell differentiation, isotype class switching, and germinal center formation.

PREPARATION & STORAGE

This monoclonal antibody preparation was purified from tissue culture supernatant via affinity chromatography. For In Vivo Ready™ (IVR) products, each preparation is also evaluated for endotoxin levels using the LAL assay. It is recommended to store the product undiluted at 4°C. Do not freeze.

APPLICATION NOTES

Tonbo Biosciences tests all of our antibodies by flow cytometry. Citations may be provided as a resource for additional applications that have not been validated by Tonbo Biosciences - please consult Materials and Methods sections for additional details about the use of any product in these publications.

REFERENCES

Clark EA, Ledbetter JA. 1986. Proc Natl Acad Sci U S A. 1986; 83(12):4494-4498. Zhang K, Clark EA, Saxon A. 1991. J Immunol. Mar 15;146(6):1836-42. Galy AH, Spits H. 1992. J Immunol. 149(3):775-782. Challa A, Pound JD, Armitage RJ, Gordon J. 1999. Allergy. Jun;54(6):576-83. Cognasse F, Chavarin P, Acquart S, Sabido O, Beniguel L, Genin C, Richard Y, Garraud O. 2005. Int J Immunopathol Pharmacol. Jan-Mar;18(1):65-74.

Tonbo Biosciences tests all antibodies by flow cytometry. Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.