

TECHNICAL DATA SHEET

Purified Mouse IgG2a Isotype Control (C1.18.4)

Catalog Number: 70-4724

PRODUCT INFORMATION

Contents: Purified Mouse IgG2a Isotype Control (C1.18.4)

Isotype: Mouse IgG2a

Concentration: 0.5 mg/mL

Clone: C1.18.4

Reactivity: N/A

Use By: 12 months from date of receipt

Storage Conditions: 2-8°C

Formulation: 10 mM NaH₂PO₄, 150 mM NaCl, 0.09% NaN₃, pH7.2

DESCRIPTION

The C1.18.4 immunoglobulin is useful as an isotype-matched control. The C.1.18.4 immunoglobulin has an unknown binding specificity and is used as an isotype control for mouse IgG2a antibodies.

PREPARATION & STORAGE

This monoclonal antibody preparation was purified from tissue culture supernatant via affinity chromatography. For In Vivo Ready™ (IVR) products, each preparation is also evaluated for endotoxin levels using the LAL assay. It is recommended to store the product undiluted at 4°C. Do not freeze.

APPLICATION NOTES

This antibody preparation has been quality-tested for flow cytometry using an appropriate cell type. The isotype control should be used at the same concentration as the primary antibody for which it is a control.

Tonbo Biosciences tests all antibodies by flow cytometry. Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.