

TECHNICAL DATA SHEET

Purified Anti-Mouse TCR beta (H57-597)

Catalog Number: 70-5961

PRODUCT INFORMATION

Contents: Purified Anti-Mouse TCR beta (H57-597)

Isotype: Armenian Hamster IgG

Concentration: 0.5 mg/mL

Clone: H57-597

Reactivity: Mouse

Formulation: 10 mM NaH₂PO₄, 150 mM NaCl, 0.09% Na₃N, pH7.2

DESCRIPTION

The H57-597 antibody is specific for the beta chain of the mouse T cell Receptor (TCR). This cell surface protein combines with a second protein chain (alpha chain) to form the alpha-beta TCR that is expressed by NK1.1+ thymocytes, NKT cells, and the majority of peripheral T cells. A small number of T cells may express an alternative heteromer of gamma and delta protein chains, known as the gamma-delta TCR. These receptors participate in a complex with CD3, and with the co-receptors CD4 or CD8, to recognize and respond to antigens bound to MHC molecules on antigen-presenting cells. Such interactions promote T cell receptor signaling (T cell activation) and can result in a number of cellular responses including proliferation, differentiation, production of cytokines or activation-induced cell death. The H57-597 antibody is used as a phenotypic marker for T cells expressing the alpha-beta TCR. It is also widely used to cross-link surface TCR and thereby mimic TCR-mediated cell activation or induction of apoptosis. The antibody does not cross-react with cells expressing the gamma-delta TCR.

PREPARATION & STORAGE

This monoclonal antibody preparation was purified from tissue culture supernatant via affinity chromatography. For In Vivo Ready™ (IVR) products, each preparation is also evaluated for endotoxin levels using the LAL assay. It is recommended to store the product undiluted at 4°C. Do not freeze.

APPLICATION NOTES

This purified format is guaranteed to be >90% pure as determined by SDS-PAGE analysis. Citations are provided as a convenience to you - please consult Materials and Methods sections for additional details about the use of any product in these publications.

REFERENCES

Berent-Maoz B, Montecino-Rodriguez E, Signer RAJ, and Dorshkind K. 2012. Blood. 119:5715-5721. (flow cytometry)Wang D, Qin H, Du W, Shen Y-W, Lee W-H, Riggs AD, and Liu C-P. 2012. Proc. Natl. Acad. Sci. 109:9493-9498. (in vitro induction of apoptosis)O'Brian RL, Taylor MA, Hartley J, Nuhsbaum T, Dugan S, Lahmers K, Aydintug MK, Wands JM, Roark CL, and Born WK. 2009. Invest. Ophthalmol. Vis. Sci. 50: 3266-3274. (immunofluorescence microscopy – OCT embedded frozen tissue)Matei IR, Gladly RA, Nutter LMJ, Carty A, Guidos CJ, and Danska JS. 2007. Blood. 109:1887-1896. (immunoprecipitation)Harada N, Shimada M, Okano S, Suehiro T, Soejima Y, Tomita Y, and Maehara Y. 2004. J. Immunol. 173:6635-6644. (in vivo T cell depletion)Kubo RT, Born W, Kappler JW, Marrack P, and Pigeon M. 1989. J. Immunol. 142: 2736-2742. (Origination of clone, immunoprecipitation, in vitro activation)

NOTE: Please choose the appropriate format for each application. Citations are provided as a convenience to you; please consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.