

TECHNICAL DATA SHEET

Purified Anti-Mouse CD279 (PD-1) (RMP1-14)

Catalog Number: 70-9982

PRODUCT INFORMATION

Contents: Purified Anti-Mouse CD279 (PD-1) (RMP1-14)

Isotype: Rat IgG2a, kappa

Concentration: 0.5 mg/mL

Clone: RMP1-14

Reactivity: Mouse

Use By: 12 months from date of receipt

Storage Conditions: 2-8°C

Formulation: 10 mM NaH₂PO₄, 150 mM NaCl, 0.09% NaN₃, pH 7.2

DESCRIPTION

The RMP1-14 antibody is specific for mouse CD279, also known as programmed death-1 (PD-1), a 55 kDa glycoprotein member of the Ig superfamily of molecules. PD-1 exists in a monomeric form that is expressed by CD4⁺ CD8⁻ thymocytes, where it participates in the processes of clonal selection, elimination of autoreactive lymphocytes, and development of tolerance. PD-1 expression is also inducible upon activation of mature T cells, where it has been proposed to interact with the co-stimulatory receptor CD80 to limit T cell activation. Two ligands for PD-1, known as PD-L1 (B7-H1) and PD-L2 (B7-DC) are differentially expressed on T and B cells, monocytes, macrophages, NK cells or dendritic cells.

PREPARATION & STORAGE

This monoclonal antibody preparation was purified from tissue culture supernatant via affinity chromatography. For In Vivo Ready™ (IVR) products, each preparation is also evaluated for endotoxin levels using the LAL assay. It is recommended to store the product undiluted at 4°C. Do not freeze.

APPLICATION NOTES

Tonbo Biosciences tests all of our antibodies by flow cytometry. Citations may be provided as a resource for additional applications that have not been validated by Tonbo Biosciences - please consult Materials and Methods sections for additional details about the use of any product in these publications.

REFERENCES

Matsumoto, K., et al. 2004. J. Immunol. 172:2530 Kanai, T., et al. 2003. J. Immunol. 171:4156.

Tonbo Biosciences tests all antibodies by flow cytometry. Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.