

TECHNICAL DATA SHEET

Goat anti-Mouse IgM (H&L) - Affinity Pure, HRP Conjugate

Catalog Number: 72-8064

PRODUCT INFORMATION

Contents: Goat anti-Mouse IgM (H&L) - Affinity Pure, HRP Conjugate

Reactivity: Mouse

Specificity: Mouse IgM (H&L)

Host: Goat

Cross Adsorption: No

Cross Reactivity: Based on IEP, no reactivity is observed to: · non-immunoglobulin mouse serum proteins

Concentration: Lyophilized

Formulation: Lyophilized, 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free

Reconstitution: Rehydrate with 1.1 ml of deionized water and let stand 30 minutes at room temperature to dissolve. Centrifuge to remove any particulates. Prepare fresh working dilution daily.

Storage Conditions: Store freeze-dried powder at 2-8 °C. Once rehydrated and in solution store at 2-8 °C.

Use By: 1 year from date of receipt. Store lyophilized material at 2-8 °C. For long term storage after reconstitution, dilute with 50% glycerol and store at -20 °C as a liquid.

Purity: > 95% based on SDS-PAGE

Purification: Affinity purified using solid phase Mouse IgM

Country of Origin: Goat serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.

PREPARATION & STORAGE

This product was Affinity purified using solid phase Mouse IgM. Store freeze-dried powder at 2-8 °C.

APPLICATION NOTES

This conjugate is suitable for Western Blot, ELISA and immunohistochemistry. The optimal working dilution should be determined by the investigator. Suggested starting dilution(κ): · 1:20 - 1:2,000 for most applications

Tonbo Biosciences tests all antibodies by flow cytometry. Citations are provided as a resource for additional applications that have not been validated by Tonbo Biosciences. Please choose the appropriate format for each application and consult Materials and Methods sections for additional details about the use of any product in these publications.

For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.