

TECHNICAL DATA SHEET

# violetFluor™ 450 Rat IgG2b Isotype Control (LTF-2)

Catalog Number: 75-4031

## PRODUCT INFORMATION

**Contents:** violetFluor™ 450 Rat IgG2b Isotype Control (LTF-2)

**Isotype:** Rat IgG2b

**Concentration:** 0.2 mg/mL

**Clone:** LTF-2

**Reactivity:** N/A

**Formulation:** 10 mM NaH<sub>2</sub>PO<sub>4</sub>, 150 mM NaCl, 0.09% NaN<sub>3</sub>,  
0.1% gelatin, pH7.2

## DESCRIPTION

The LTF-2 immunoglobulin is useful as an isotype-matched control. The LTF-2 immunoglobulin has an unknown binding specificity and is used as an isotype control for rat IgG2b antibodies.

## PREPARATION & STORAGE

This monoclonal antibody preparation was purified from tissue culture supernatant via affinity chromatography. For In Vivo Ready™ (IVR) products, each preparation is also evaluated for endotoxin levels using the LAL assay. It is recommended to store the product undiluted at 4°C. Do not freeze.

## APPLICATION NOTES

This antibody preparation has been quality-tested for flow cytometry using an appropriate cell type. The isotype control should be used at the same concentration as the primary antibody for which it is a control.

violetFluor™ 450 dye is excited by the violet (405 nm) laser and has a peak emission of 450 nm. The most common band pass filters for this dye are 440/40 or 450/50. violetFluor™ 450 can be used as an alternative for Pacific Blue®, BD Horizon™ V450 or eFluor® 450.

## For Research Use Only.

Not for use in diagnostic or therapeutic procedures. Not for resale. Not for distribution without written consent. Tonbo Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Tonbo Biosciences, Tonbo Biosciences Logo and all other trademarks are the property of Tonbo Biotechnologies Corporation. © 2013 Tonbo Biosciences.