



## TECHNICAL DATA SHEET

### Shutdown Fluid (1X)

Catalog Number: TNB-4601-L1500

## PRODUCT INFORMATION

**Contents:** Shutdown Fluid (1X)

**Use By:** 6 months from date of receipt

**Storage Conditions:** 18-25°C

**Formulation:** 1X formulation containing 2-phenoxyethanol

## DESCRIPTION

Cytek Biosciences Shutdown Fluid (1X) is a pre-filtered 2-phenoxyethanol solution that can be used in the shutdown cycle of Cytek Biosciences cytometers to maintain instrument cleanliness during long-term storage. Each bottle provides 1.5 liter of ready-to-use Shutdown Fluid.

## USE

Shutdown Fluid is provided as a ready-to-use 1X solution for long-term storage of Cytek Biosciences analyzers. Add the contents of the Shutdown Fluid (1X) container to a clean auxiliary sheath tank and run a Long Clean cycle.

## STORAGE

Product should be stored at 18-25°C. Do not freeze.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale or distribution without Cytek's written consent. Cytek and the Cytek logo are trademarks of Cytek Biosciences, Inc. © 2022 Cytek Biosciences, Inc. All rights reserved.